

Chemistry Department

Standard Operating Procedure

Title: Pyrophoric Chemicals

Pyrophoric chemicals are liquids and solids that will ignite spontaneously in air at about 130 F. Titanium dichloride and phosphorus are example of pyrophoric solids; tributylaluminum and related compounds are examples of pyrophoric liquids.

Decontamination procedures

Personnel: Wash hands and arms with soap and water immediately following any skin contact with pyrophoric chemicals.

Emergency procedure

Emergency procedures that address response actions to fires, explosions, spills, injury to staff, or the development of sign and symptom of overexposure must be developed. The procedures should address as a minimum the following:

Who to contact: (University police, and Office of Environmental Health and Safety, Principal investigator of the laboratory including evening phone number)

The location of all safety equipment (showers, eye wash, fire extinguishers, etc.)

The method used to alert personnel in nearby areas of potential hazards

Special spill control materials required by the type of pyrophoric chemicals handled in the laboratory

Eye protection

Eye protection in the form of safety glasses must be worn at all times when handling pyrophoric chemicals. Ordinary (street) prescription glasses do not provide adequate protection. (Contrary to popular opinion these glasses cannot pass the rigorous test for industrial safety glasses.) Adequate safety glasses must meet the requirements of the Practice for Occupational and Educational Eye and Face Protection (ANSI Z.87. 1 1989)

Developed on
7/18/2002

Developed by
Gary Johnson
July, 2002

Chemistry Department

Standard Operating Procedure

and must be equipped with side shields. Safety glasses with side shields do not provide adequate protection from splashes; therefore, when the potential for splash hazard exists other eye protection and/or face protection must be worn.

Eyewash

Where the eyes or body of any person may be exposed to pyrophoric chemicals, suitable facilities for quick drenching or flushing of the eyes and body shall be provided within the work area for immediate emergency use. Bottle type eyewash stations are not acceptable.

Fume hood

Many pyrophoric chemicals release noxious or flammable gases and should be handled in a hood. In addition some pyrophoric materials are stored under kerosene (or other flammable solvents), therefore the use of a fume hood is required to prevent the release of flammable vapors in the laboratory. Glove boxes may be also be used.

Glove (dry) box

Glove boxes may be used to handle pyrophoric chemicals if inert or dry atmospheres are required.

Gloves

Gloves should be worn when handling pyrophoric chemicals. Disposable latex or nitrile gloves provide adequate protection against accidental hand contact with small quantities of most laboratory chemicals. Lab workers should contact the department store room for advice on chemical resistant glove selection when direct or prolonged contact with hazardous chemicals is anticipated.

Chemistry Department

Standard Operating Procedure

Hazard assessment

Hazard assessment for work involving pyrophoric chemicals should thoroughly address the issue of fire safety (including the need for Class D fire extinguishers), proper use and handling techniques, chemical toxicity, storage, and spill response.

Protective apparel

Lab coats, closed toed shoes and long sleeved clothing should be worn when handling pyrophoric chemicals. Additional protective clothing should be worn if the possibility of skin contact is likely.

Safety shielding

Safety shielding is required any time there is a risk of explosion, splash hazard or a highly exothermic reaction. All manipulations of pyrophoric chemicals that pose this risk should occur in a fume hood with the sash in the lowest feasible position. Portable shields, which provide protection to all laboratory occupants, are acceptable.

Safety shower

A safety or drench shower should be available in a nearby location where the pyrophoric chemicals are used.

Signs and labels

Containers: All pyrophoric chemicals must be clearly labeled with the correct chemical name. Handwritten labels are acceptable; chemical formulas and structural formulas are not acceptable.

Special storage

Pyrophoric chemicals should be stored under an atmosphere of inert gas or under kerosene as appropriate. Do not store pyrophoric chemicals with flammable materials or in

Developed on
7/18/2002

Developed by
Gary Johnson
July, 2002

Chemistry Department

Standard Operating Procedure

a flammable-liquids storage cabinet. Store these materials away from sources of ignition. Minimize the quantities of pyrophoric chemicals stored in the laboratory.

Never return excess chemicals to the original container. Small amounts of impurities introduced into the container may cause a fire or explosion.

Special ventilation

Always attempt to handle pyrophoric chemicals in a fume hood or glove box. If your research does not permit the handling of pyrophoric chemicals in a fume hood or glove box you must consult with the Principle Investigator to review the adequacy of any special ventilation decisions.

Spill response

Anticipate spills by having the appropriate clean up equipment on hand. The appropriate clean up supplies can be determined by consulting the material safety data sheet. This should occur prior to the use of any pyrophoric chemicals. Spill control materials for pyrophoric chemicals are designed to be inert and will not react with the reagent.

In the event of a spill alert personnel in the area that a spill has occurred. Do not attempt to handle a large spill of pyrophoric chemicals. Turn off all ignition sources and vacate the laboratory immediately. Call for assistance.

- Office of Environmental Health & Safety 335-3041 or 911
- University Police 911
- This is a 24-hour service.

Remain on the scene, but at a safe distance, to receive and direct safety personnel when they arrive.

Vacuum protection

Evacuated glassware can implode and eject flying glass, and splattered chemicals. Vacuum work involving pyrophoric chemicals must be conducted in a fume hood or isolated in an acceptable manner.

Developed on
7/18/2002

Developed by
Gary Johnson
July, 2002

Chemistry Department

Standard Operating Procedure

Mechanical vacuum pumps must be protected using cold traps and, where appropriate, filtered to prevent particulate release. The exhaust for the pumps must be vented into an exhaust hood. Vacuum pumps should be rated for use with pyrophoric chemicals.

Waste disposal

All materials contaminated with pyrophoric chemicals should be disposed of as hazardous waste. Contact the chemistry department storeroom to dispose of wastes. See the Waste Disposal SOP for collection, storing, submitting and labeling instructions. These wastes may pose a flammability risk and should not remain in the laboratory for any extended time.